

## BL21(DE3) Chemically Competent Cell

Cat. No. CD601

**Storage:** at -70°C for six months. Do not store in liquid nitrogen.

### Description

BL21(DE3) Chemically Competent Cell is specifically designed for chemical transformation of DNA. It permits a transformation efficiency of over  $10^7$  cfu/ $\mu$ g DNA (tested by pUC19 plasmid DNA).

### Genotype

F<sup>-</sup> *ompT hsdS*(r<sub>B</sub><sup>-</sup>m<sub>B</sub><sup>-</sup>) *gal dcm*(DE3)

### Features

- Transformation efficiency:  $>10^7$  cfu/ $\mu$ g (pUC19 DNA).
- DE3 strains contains the  $\lambda$ DE3 lysogen that carries the gene for T7 RNA polymerase.
- Suitable for T7 and T7lac such as pET, *pEASY*<sup>®</sup>.
- Suitable for high expression of non-toxic protein.
- Control plasmid I (Amp<sup>r</sup>) is used for detection of expression function of cell. The protein size is about 25 kDa.

### Procedures

- Equilibrate a water bath to 42°C.
- Warm a vial of SOC medium or LB medium to room temperature. Warm selective plates at 37°C for 30 minutes.
- Thaw a vial of 100  $\mu$ l of BL21(DE3) Chemically Competent Cell on ice, aliquot 50  $\mu$ l of the cells into a prechilled 1.5 ml tube, add target DNA (1 to 5  $\mu$ l) into the tube. Do not mix by pipetting up and down. Incubate the cells on ice for 30 minutes.
- Heat-shock the cells for 45 seconds at 42°C without shaking. Immediately transfer the tube to ice. Incubate on ice for 2 minutes without shaking.
- Add 500  $\mu$ l of prewarmed SOC medium or LB medium (without antibiotic) into the tube, mix well and shake at 37°C for 1 hour at 200 rpm.
- Spread 20 to 200  $\mu$ l from each transformation vial on a prewarmed selective plate. The remaining can be stored at 2-8°C and plated the next day if needed.
- Invert the plate and incubate at 37°C overnight.
- Select colonies and analyze by restriction enzyme digestion, PCR, or sequencing.

### Notes

- Higher efficiency transformation can be achieved by transforming cells immediately following thawing.
- Avoid repeated thawing.
- Gentle handling is required for the entire procedure.

FOR RESEARCH USE ONLY