



42<sup>°C</sup> 65<sup>°C</sup>

## TransScript<sup>®</sup>-Uni One-Step gDNA Removal and cDNA Synthesis SuperMix

TransScript<sup>®</sup>-Uni RT is designed by genetic modification. It provides a broad range of reaction temperature (42°C-65°C) and ultra-high thermostability. This supermix achieves simultaneous genomic DNA removal and cDNA synthesis in one tube. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds. The suggested reaction temperature is 50°C.

- Ultra-high thermostability: reaction temperature at 42°C-65°C.
- Simultaneous genomic DNA removal and cDNA synthesis in one tube to minimize RNA contamination.
- After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated simultaneously. Compared with traditional treatment with DNase I before reverse transcription, inactivation of DNase I is not required and degradation of RNA is avoided.
- cDNA fragment up to 20 kb.

We provide limited quantity of free samples, please contact with us for ordering.

Product Name	Cat. No.	Specification
TransScript® -Uni One-Step gDNA Removal and cDNA Synthesis SuperMix	AU311-02	50 rxns×20 µl System
	AU311-03	100 rxns×20 µl System

### Tips

TransScript® -Uni One-Step gDNA Removal and cDNA Synthesis SuperMix is stored at -20°C for one year.

### Notes

Further details can be found in the customer contact below. For any questions, please contact with us by WeChat or call Customer Service at +86-10-57815027/57815087 or contact with our distributors.



## About us

- ▶ We provide high quality products for life science.
- ▶ We provide simple and efficient methods.
- ▶ We provide new concept for scientific research.

## TRANSGEN BIOTECH CO., LTD.




Website [www.transgenbiotech.com](http://www.transgenbiotech.com)

E-mail [info@transgenbiotech.com](mailto:info@transgenbiotech.com)

Customer Service +86-400-898-0321

Phone +86-10-57815027

B to B Website <http://transgen.en.alibaba.com>

 Facebook
  Twitter
  LinkedIn
  Skype



WeChat