



*The **BEST** for
Life Science*

6×Protein Loading Buffer

Cat. No. DL101

Storage: at -20°C for two years

Description

6×Protein Loading Buffer is especially formulated for protein sample preparation to be used in SDS-PAGE. Prior to loading, add appropriate volume of 6×Protein Loading Buffer to protein sample to make its working concentration at 1×.

FOR RESEARCH USE ONLY

Website www.transgenbiotech.com

E-mail info@transgenbiotech.com

Customer Service +86-400-898-0321

Phone +86-10-57815027